

Contents

1. Description

1.1 Background information

1.2 Applications

2. Protocol

2.1 Preparation

2.2 Instructions for use

1. Description

Components	100 mL MACS® Tissue Storage Solution
Sterility	Sterile filtered.
Storage	Store protected from light at 2–8 °C. Do not freeze. The expiration date is indicated on the label.

1.1 Background information

The storage of organs or tissue samples may be required if material has to be collected from multiple origins, at different time points, or if it cannot be directly processed. During this storage period necrosis and apoptosis must be avoided and the activation and pluripotency status of cells has to be preserved. The MACS® Tissue Storage Solution has been developed to allow an optimized storage of fresh organ and tissue samples for at least 48 hours without showing background effects like cell activation or apoptosis induction. It has been tested and validated on a variety of human and murine tissues including tumors, skin, heart, spleen, brain, and skeletal muscle.

1.2 Applications

- Storage of fresh organ and tissue samples
- Shipping of fresh organ and tissue samples

2. Protocol

2.1 Preparation

- ▲ The MACS Tissue Storage Solution is formulated without antibiotic or antifungal components. If required, supplement with these ingredients at standard cell culture concentrations prior to use.
- ▲ All procedures should be carried out using sterile techniques in a biosafety cabinet.
- ▲ Pre-chill MACS Tissue Storage Solution to 2–8 °C before use.

2.2 Instructions for use

1. Transfer sample directly into the pre-chilled MACS Tissue Storage Solution (2–8 °C). If the sample has already been placed in other buffers or media, replace the complete liquid with the pre-chilled MACS Tissue Storage Solution. The sample must be completely covered with the MACS Tissue Storage Solution.
2. (Optional) For shipping purpose, transfer sample into the MACS Tissue Storage Solution in a closable tube. Avoid air bubbles to prevent air contact for prolonged time periods.
3. Store or ship sample at 2–8 °C.
4. At the end of storage, remove sample from the MACS Tissue Storage Solution and process directly with downstream application, e.g., dissociation into single-cell suspensions by using a gentleMACS Dissociator and a MACS Dissociation Kit, e.g., gentleMACS Octo Dissociator with Heaters (# 130-096-427) and Tumor Dissociation Kit, mouse (# 130-096-730).

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

gentleMACS, MACS, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.