

Contents

- 1. Description
 - 1.1 Background information
 - 1.2 Applications
 - 1.3 Reagent and instrument requirements
- 2. Protocol
 - 2.1 Sample preparation
 - 2.2 T cell activation and expansion
 - 2.3 Immunofluorescent staining
- 3. Examples of T cell activation and expansion using the T Cell TransAct

1. Description

This product is for research use only.

T Cell TransAct human Components

Components	i Gen Hunshet, numun			
	Product	Order no.		
	1×2 mL T Cell TransAct, human	130-128-758		
	2×2 mL T Cell TransAct, human	130-111-160		
Capacity	1×2 mL T Cell TransAct, human is sufficient to activate and expand up to 2×10^8 enriched T cells or up to 4×10^8 peripheral blood mononuclear cells (PBMCs) or , when used at recommended titer of 1:100. 2×2 mL T Cell TransAct, human is sufficient to activate and expand up to 4×10^8 enriched T cells or up to 8×10^8 peripheral blood mononuclear cells (PBMCs), when used at recommended titage of 1:100			
Product format	roduct format Polymeric nanomatrix c humanized CD3 and CD28 a in phosphate-buffered-saline (F			
Storage	Store protected from light at $2-8$ °C. Do not freeze. The expiration date is indicated on the vial label.			

1.1 Background information

The T Cell TransAct has been designed to activate and expand enriched T cell populations or human resting T cells from periphal blood mononuclear cells (PBMCs). T cell expansion is achieved by culturing for up to 14 days. For longer cultivation restimulation after 14 days is necessary.

Polyclonal T cell expansion can be used when increased numbers of T cells are required or when T cells are activated to enable gene modification.

Due to the nanomatrix of the T Cell TransAct, it can be sterile filtered and excess reagent can be removed by simple replacement of supernatant or by a washing step, e.g., centrifugation.

T Cell TransAct™

The recommended titers have been found to efficiently stimulate the majority of T cell subsets, however, for special applications it is recommended to experimentally determine the optimal stimulation titer. Over-activation of T cells carries a risk of activation-induced cell death.

The T Cell TransAct has been developed in combination with the TexMACS[™] Medium and Human IL-2 IS or Human IL-7 and Human IL-15.

1.2 Applications

human

The T Cell TransAct is intended for the in vitro stimulation and expansion of purified T cell populations of, for example, untouched T cells isolated with the Pan T Cell Isolation Kit, human, as well as of human T cells from hematological cell populations (e.g. PBMCs).

1.3 Reagent and instrument requirements

- TexMACS™ Medium, (# 130-097research grade 196) supplemented with Human IL-2 IS, premium grade (#130-097-744) or Human IL-7, premium grade (# 130-095-361) and Human IL-15, premium grade (# 130-095-762).
- Buffer for flow cytometric analysis: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS® BSA Stock Solution (#130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (2-8 °C).

▲ Note: EDTA can be replaced by other supplements such as anticoagulant citrate dextrose formula-A (ACD-A) or citrate phosphate dextrose (CPD). Buffers or media containing $\mathrm{Ca}^{^{2+}}$ or $\mathrm{Mg}^{^{2+}}$ are not recommended for use.

- Fluorochrome-conjugated antibodies for flow cytometric analysis, e.g., CD4 Antibody, anti-human, VioBlue®, REAfinity[™], CD8 Antibody, anti-human, VioGreen[™], REAfinity, CD25 Antibody, anti-human, PE, REAfinity, and CD69 Antibody, anti-human, APC, REAfinity. For more information about fluorochrome-conjugated antibodies refer to www.miltenyibiotec.com/antibodies.
- (Optional) Pan T Cell Isolation Kit, human (# 130-096-535)
- (Optional) Propidium Iodide Solution (# 130-093-233) or 7-AAD Staining Solution (# 130-111-568) for flow cytometric exclusion of dead and apoptotic cells.

2. Protocol

▲ All steps in the protocol have to be performed under steril conditions.

▲ Excess of T Cell TransAct is removed by simple replacement of supernatant or by a washing step, e.g., centrifugation (at least 10-fold reduction) 2-3 days after initial stimulation. Performing a

www.miltenyibiotec.com

Miltenvi Biotec B.V. & Co. KG

Friedrich-Ebert-Straße 68, 51429 Bergisch Gladbach, Germany Phone +49 2204 8306-0, Fax +49 2204 85197

^{140-004-952.05} macsde@miltenyi.com

washing step earlier may result in reduced T cell proliferation.

▲ Activated T cells can be transduced 1–2 days after activation. The optimal virus titer has to be defined before and depends on the viral vector used. The T Cell TransAct can be used in combination with retro- or lenti-viral transduction.

▲ Presence of residual EDTA, e.g., when using medium containing EDTA for T cell purification, will hamper T cell stimulation. Ensure extensive removal of EDTA (i.e. over 200-fold reduction) prior to T cell stimulation with the T Cell TransAct.

2.1 Sample preparation

When working with anticoagulated peripheral blood or buffy coat, peripheral blood mononuclear cells (PBMCs) should be isolated by density gradient centrifugation, for example, using Ficoll-Paque[™].

▲ Note: To remove platelets after density gradient separation, resuspend cell pellet in buffer and centrifuge at 200×g for 10–15 minutes at 20 °C. Carefully aspirate supernatant. Repeat washing step.

For details refer to the protocols section at www.miltenyibiotec.com/ protocols.

For the isolation of purified T cells use, for example, the Pan T Cell Isolation Kit, human.

2.2 T cell activation and expansion

The protocol has been optimized for gentle and efficient activation and expansion of purified T cells and PBMCs by using a titer of 1:100.

Purified T cells should be activated at an optimal surface density of 1×10^6 cells per cm² (table 1) and PBMCs with up to 2×10^6 per cm².

Culture plate	Growth area	Max. working volume	Total T cell number	T Cell TransAct to add per well
96 well	0.31 cm ²	0.2 mL	0.3×10 ⁶	2 µL
48 well	1 cm ²	1 mL	1×10 ⁶	10 µL
24 well	2 cm ²	2 mL	2×10 ⁶	20 µL
12 well	4 cm ²	4 mL	4×10 ⁶	40 µL
6 well	10 cm ²	5 mL	5×10 ⁶	50 μL

 Table 1: Optimal surface densitiy when working with purified T cells.

Volumes given below are for the stimulation in a 48-well plate of up to 1×10^6 purified T cells or up to 2×10^6 PBMCs in a total volume of 990 µL TexMACS Medium supplemented with 20 IU/mL Human IL-2 or 155 U/mL Human IL-7 and 290 U/mL Human IL-15. When working with fewer than 10⁶ cells, use the same volumes as indicated in table 1. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly.

Activation in a 48-well plate

- 1. Determine cell number.
- 2. Resuspend cells in 990 µL supplemented TexMACS Medium.
- 3. Add 10 µL of the T Cell TransAct.
- 4. Incubate at 37 °C, 5% CO₂ for up to 3 days.
 ▲ Note: Inspect culture daily, and add fresh medium if required.
- 5. Remove residual reagent 2–3 days after initial activation by either replacing 900 μ l of supernatant with fresh supplemented TexMACS Medium or by centrifugation at 300×g for

10 minutes and aspirate supernatant completely.

6. Add 1 mL fresh supplemented TexMACS Medium and incubate at 37 °C, 5% CO₂.

Expansion

- 1. Split cell suspension every 2 days into two equal parts and add fresh supplemented TexMACS Medium.
- Incubate at 37° C, 5% CO₂.
 ▲ Note: For optimal expansion of T cells a daily inspection of culture is required. It might be necessary to split culture more or less frequently than every day.
- 3. At day 14 proceed to downstream application, e.g., analysis of cells.
- 4. (Optional) T cells can be further expanded by reapplying T Cell TransAct to the culture. However, for restimulation it is recommended to use a titer of 1:500.

2.3 Immunofluorescent staining

▲ Work fast, keep cells cold, and use pre-cooled solutions. This will prevent capping of antibodies on the cell surface and non-specific cell labeling.

▲ It is recommended to stain 10^6 cells per sample. When working with up to 10^7 cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10^7 nucleated cells, use twice the volume of all indicated reagent volumes and total volumes).

▲ The recommended incubation temperature is 2–8 °C. Higher temperatures and/or longer incubation times may lead to non-specific cell labeling. Working on ice may require increased incubation times.

▲ Upon stimulation, expression of CD3 will be transiently downregulated. Thus, the staining of CD3 on the cell surface of activated cells might be affected.

- 1. Determine cell number.
- Wash cells by adding 1-2 mL of buffer per 10⁶ cells and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
- 3. Add each staining antibody, e.g., CD4-VioBlue, CD8-VioGreen, CD25-PE, and CD69-APC according to manufacturer's recommendations.
- 4. Mix well and incubate for 10 minutes in the dark in the refrigerator (2–8 °C).
- Wash cells by adding 1-2 mL of buffer per 10⁶ cells and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
- 6. Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry or fluorescence microscopy.

Unless otherwise specifically indicated, all Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use.

3. Examples of T cell activation and expansion using the T Cell TransAct

A) Example of a T cell activation

Human purified T cells were isolated using the Pan T Cell Isolation Kit and activated for 48 hours using the T Cell TransAct (titer 1:100) in TexMACS Medium supplemented with Human IL-2 (20 IU/mL). The negative control experiment was performed without adding the T Cell TransAct. Cells were fluorescently stained using CD25-PE and CD69-APC and analyzed by flow cytometry using the MACSQuant[®] Analyzer. CD4-VioBlue was used for selection of T helper cells and CD8-VioGreen was used for selection of cytotoxic T cells. Dead cells and debris were excluded from the analysis based on scatter signals and propidium iodide fluorescence.

Negative control





B) Expansion of pan T cells after activation

Isolated Pan T cells were labeled with carboxyfluorescein succinimidyl ester (CFSE) and stimulated with the T Cell TransAct. T cells were cultured at a density of 1×10^6 cells per cm² in supplemented TexMACS Medium supplemented with Human IL-2 (20 IU/mL). Proliferation analysis was done by flow cytometry via the detection of the CFSE dilution 7 days after stimulation. Non-stimulated pan T cells act as negative control.



Refer to **www.miltenyibiotec.com** for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shell file stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, MACS, MACSQuant, the Miltenyi Biotec logo, REAfinity, TexMACS, TransAct, VioBlue, and VioGreen are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Ficoll-Paque is a trademark of GE Healthcare companies

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.

Unless otherwise specifically indicated, all Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use.